



chem-supply

Buffer Solutions



Can't afford calibration mistakes?
 Tired of erroneous pH readings?
 Using buffers that aren't up to the job?



WE HAVE THE SOLUTIONS.

Scharlau offers buffer solutions that are **precise, accurate, and reliable.**

You can use Scharlau buffers with **complete confidence** knowing they are:

- **NIST** traceable
- certificated to an uncertainty of just **±0.01 pH** unit
- available in **500mL** and **1L** pack sizes
- available **colour-coded** (for easy identification), and also as 'clear' (colourless)
- supplied *individually wrapped* together with its own **Certificate of Analysis**
- supplied with a **temperature adjustment** scale on each bottle (& CofA)
- have a specified **2 year shelf life** so you know your product has not expired
- batch values are determined by a **five-point calibration** measurement

Product Code	Description
SO2004	Buffer solution pH 4.00 (20 °C), red
SO2007	Buffer solution pH 7.00 (20 °C), yellow
SO2010	Buffer solution pH 10.00 (20 °C), blue
SO1004	Buffer solution Ph 4.00 (20 °C) clear
SO1007	Buffer solution pH 7.00 (20 °C) clear
SO1010	Buffer solution pH 10.00 (20 °C) clear

Contact your local Chem-Supply Distributor or Chem-Supply for more details.



38-50 Bedford Street, Gillman SA 5013
 P: + 61 (0) 8 8440 2000
 F: + 61 (0) 8 8440 2001
www.chemsupply.com.au