

## Dichloromethane, Pesticide

CH <sub>2</sub> Cl <sub>2</sub>	FW. 84.93	Density 1 L = 1.330 Kg.
CAS-No.	75-09-2	Melting Point - 95 °C
Code	PC1040A	Boiling Point 40 °C

## Specifications

Assay (by GC.)	99.9% min.
Identity (IR)	Passes test
Color (APHA)	10 max.
Water (by Coulometry)	0.01% max.
Acidity (mEq./g.)	0.0003 max.
Alkalinity (mEq./g.)	0.0002 max.
Residue on Evaporation	0.0003% max.
ECD (as lindane standard) Single impurity peak	10 ng/L

Stabilized with about 50 ppm amylene.