

Methanol, Purge and Trap

CH₃OH FW. 32.04
CAS-No. 67-56-1
Code [PT1115](#)

Density 1 L = 0.790 Kg.
Melting Point - 98 °C
Boiling Point 64.5 °C

Specifications

Assay (by GC.)	99.9%	min.
Identity (IR)	Passes test	
Color (APHA)	10	max.
Water (by Coulometry)	0.05%	max.
Acidity (mEq./g.)	0.0005	max.
Alkalinity (mEq./g.)	0.0002	max.
Residue on Evaporation	0.0002%	max.
Methyl ethyl ketone (by GC/MS)	50	ppb max.
Toluene (by GC/MS)	50	ppb max.

Product passed through 0.2 micron final filter.